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
Indications

CombiHaler® Solo is an inhalation chamber intended for invasive and non invasive mechanical ventilation circuits. This device was designed for the administration of aerosol drugs for mechanically ventilated adult patients. CombiHaler® Solo is intended for use with a pressurized metered dose Inhaler (MDI). CombiHaler® Solo does not weight down on the ventilation circuit and allows monitoring of medication administration.

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COMPOSITION

- Cone connector 22 mm
- Socket connector 22 mm
- One opening for MDI
- One sealed cap
- Nozzle


Prerequisites before use

- The bag must be opened only when you intend to use CombiHaler® Solo. Once out of its packaging, this product must be used immediately.
- When starting, it is necessary to ensure that :
 - The device does not contain any foreign objects, to avoid any risk of damage to the patient's airways or obstruction of the medical device.
 - The MDI nozzle correctly ensures its function (visually test with a puff of drugs), in order to avoid a risk of incorrect delivery of the medicine.
- Outside periods of drug administration, it must be ensured that the cap is properly closed, to guarantee the tightness of the device and to avoid any risk of leakage.

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


Instructions for use

Checking the operation of the nozzle

Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister.

Operate the MDI and inspect the nozzle for proper function. Look at the medication jet.

Insertion in the circuit


In a traditional adult dual-limb mechanical ventilation circuit, insert CombiHaler® Solo on the inspiratory branch of the circuit. The cone connector (a) must be plugged to the inspiratory branch. The socket connector (b) must be plugged to the Y piece. The arrow on the device must point towards the patient.

Drug Administration by pressurized dose metered

Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Open the cap from the opening (c) and insert the canister.

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Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI.

Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c).

Bénéfices cliniques – performances du dispositif / Clinical benefits – device performance

Performance *in vitro* / In Vitro performance

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The CombiHaler® Solo inhalation chamber was evaluated in an adult mechanical ventilation circuit with a ventilator (Dräger, Evita 2 Dura) set to volume-controlled mode. 10 doses of 100 μg salbutamol (Ventolin®, GSK) were activated in the device during the inspiratory phase at intervals of 1 minute. An endotracheal tube (ET) with an internal diameter of 7.5 mm was connected to the Y-piece and to a test lung model (SmartLung Adult, IMT Medical). The deposited dose was recovered on a filter positioned between the ET and the test lung and was quantified by UV spectrophotometry.

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Storage – handling – disposal

- The recommended lifespan of CombiHaler® Solo is equivalent to the recommended service life for the mechanical ventilation circuits to which it is associated, without exceeding 7 days, or sooner if a visual examination indicates that it is necessary to make its change.
- CombiHaler® Solo must be disposed of at the same time as the mechanical ventilation circuit, in accordance with the instructions of the facility.
- CombiHaler® Solo must not be used on more than one circuit and must be destroyed under the same conditions as the circuit on which it was mounted.
- Any total or partial replacement of the circuit must result in the replacement of CombiHaler® Solo.
- Store away from moisture.

Precautions for use

- CombiHaler® Solo is a single patient use product. Do not clean and sterilize. Do not reuse to avoid the risk of contamination between patients.
- To ensure that the device remains in place, the CombiHaler® Solo must not be dismantled. If the CombiHaler® Solo is dismantled, it must not be reused and must be replaced.

Warnings – contraindications

- This device requires to be used by qualified staff to take care of patients who are under respiratory assistance.
- CombiHaler® Solo is designed to be inserted into an adult mechanical ventilation system.
- Any heat and moisture exchanger, located between the spacer and the patient, must be temporarily removed during the use of CombiHaler® Solo.
- Do not use CombiHaler® Solo with a pressure over 35 mbar without careful and permanent attention.

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The trademark CombiHaler® belongs to OptimHal-ProtécSom SAS.
OptimHal-ProtécSom SAS,
CombiHaler® est une marque de commerce d'OptimHal-ProtécSom SAS.
Brevets déposés conjointement par l'Université de Tours, l'INSERM et OptimHal-ProtécSom SAS.
Patents belongs to Tours University, INSERM, and OptimHal-ProtécSom SAS together.