# CombiHaler®

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### COMPOSITION

- Cone connector 22 mm
- Socket connector 22 mm
- One opening for MDI
- One opening for vibrating mesh nebulizer
- Two sealed caps
- Nozzle
- A Respimat® inhaler adapter can also be supplied

### Indications

CombiHaler is an inhalation chamber intended for invasive and non invasive mechanical ventilation circuits. This device was designed for the administration of aerosol drugs for mechanically ventilated patients. It allows the use of both metered dose inhaler and vibrating mesh nebulizer. CombiHaler® does weight down only a little on the ventilation circuit and allows monitoring of medication administration.


### Prerequisites before use

- The bag must be opened only when you intend to use CombiHaler®. Once out of its packaging, this product must be used immediately.
- When starting, it is necessary to ensure that :
  - The device does not contain any foreign objects, to avoid any risk of damage to the patient's airways or obstruction of the medical device.
  - The MDI nozzle correctly ensures its function (visually test with a puff of drugs), in order to avoid a risk of incorrect delivery of the medicine.
- Outside periods of drug administration, it must be ensured that the caps are properly closed, to guarantee the tightness of the device and to avoid any risk of leakage.

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




3.

## Instructions for use

### Checking the operation of the nozzle

Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister.

Operate the MDI and inspect the nozzle for proper function. Look at the medication jet.

### Insertion in the circuit

In a traditional adult dual-limb mechanical ventilation circuit, insert CombiHaler® on the inspiratory branch of the circuit. For example, for the device to be inserted just before the Y piece, the cone connector (a) must be plugged to the inspiratory branch. The socket connector (b) must be plugged to the Y piece. The arrow on the device must point towards the patient.

### Drug Administration by pressurized dose metered

Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Open the cap from the opening (c) and insert the canister.


Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI.

Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c).

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## Administration of drugs by vibrating mesh nebulizer


Remove the cap from the orifice (d) and insert the vibrating mesh nebulizer. To avoid an opening of the circuit this step could be done before insertion of the spacer in the circuit. Deposit the drug solution into the nebulizer reservoir. For proper operation of the nebulizer, follow the instructions provided by the manufacturer.

The vibrating mesh nebulizers Aeroneb Solo® (Aerogen) can remain on the circuit during its service life as recommended by the manufacturer.

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## Drug administration with the RespiMat® inhaler

Remove the cap from the orifice (d) and insert the adapter RespiMat®. Insert the RespiMat inhaler into the adapter. To operate the inhaler, follow the instructions provided by the manufacturer.

Repeat step 1 according to the number of doses required. Follow the instructions of the manufacturer of the MDI for the waiting time between each dose. As soon as the medication is complete, remove the inhaler and adapter and put the cap back on the orifice (d).

## Performance *in vitro* / In Vitro performance

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The CombiHaler® inhalation chamber was evaluated in an adult mechanical ventilation circuit with a ventilator (Dräger, Evita 2 Dura) set to volume-controlled mode. 10 doses of 100 µg salbutamol (Ventolin® , GSK) were activated in the device during the inspiratory phase at intervals of 1 minute. An endotracheal tube (ET) with an internal diameter of 7.5 mm was connected to the Y-piece and to a test lung model (SmartLung Adult, IMT Medical). The deposited dose was recovered on a filter positioned between the ET and the test lung and was quantified by UV spectrophotometry.

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## Storage – handling – disposal

- The recommended lifespan of CombiHaler® is equivalent to the recommended service life for the mechanical ventilation circuits to which it is associated, without exceeding 7 days, or less if a visual examination indicates that it is necessary to make its change.
- CombiHaler® must be disposed of at the same time as the mechanical ventilation circuit, in accordance with the instructions of the facility.
- CombiHaler® must not be used on more than one circuit and must be destroyed under the same conditions as the circuit on which it was mounted.
- Any total or partial replacement of the circuit must result in the replacement of CombiHaler®.
- Store away from moisture.

## Precautions for use

- CombiHaler® is a single patient use product. Do not clean and sterilize. Do not reuse with others patients to avoid the risk of contamination.
- To ensure that the device remains in place, the CombiHaler® must not be dismantled. If the CombiHaler® is dismantled, it must not be reused and must be replaced.

## Warnings – contraindications

- This device requires to be used by qualified staff to take care of patients who are under respiratory assistance.
- CombiHaler® is designed to be inserted into an adult mechanical ventilation system.
- Any heat and moisture exchanger, located between the spacer and the patient, must be temporarily removed during aerosol drug delivery.
- Do not use CombiHaler® with a pressure over 35 mbar without careful and permanent attention.

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The trademark CombiHaler® belongs to OptimHal-ProtécSom SAS.
**OptimHal-ProtécSom SAS,**
CombiHaler® est une marque de commerce d'OptimHal-ProtécSom SAS.
**Brevets déposés conjointement par l'Université de Tours, l'INSERM et OptimHal-ProtécSom SAS.**
**Patents belongs to Tours University, INSERM, and OptimHal-ProtécSom SAS together.**