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

Low volume (60 mL) inhalation chamber for mechanical ventilation circuit

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



MinimHal® is an inhalation chamber intended for adult or pediatric patients with invasive or non-invasive mechanical ventilation who have been prescribed aerosolized treatment from a pressurized metered-dose inhaler and/or a vibrating mesh nebulizer used simultaneously or not.

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- MinimHal® should be used by clinicians qualified in the care of mechanically ventilated patients.
- Please read the instructions for use of MinimHal®, the metered dose inhaler and the nebulizer and keep them available at all times of use.

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- a. Male connector 15 mm
- b. Female connector 22 mm
- c. One opening for MDI
- d. One opening for vibrating mesh nebulizer
- e. Two sealed caps
- f. One nozzle
- g. A 15 mm female to 22 mm male converter that allows the use of MinimHal® on a pediatric circuit or on an adult circuit.
- h. A RespiMat® inhaler adapter can also be supplied.

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When putting the device into operation, it is imperative to ensure that:

- The device does not contain any foreign body, to avoid any risk of damage to the patient's airways or obstruction of the medical device.
- The caps are correctly positioned when MDI and/or nebulizer is not used.

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Checking the operation of the nozzle
Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister (Fig.2).
Operate the MDI and check that the nozzle correctly ensures its function (visually test with a puff of drugs) (Fig.3).

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For insertion of MinimHal® on the inspiratory limb of a pediatric dual limb circuit:
The 15mm male connector must be connected to the inspiratory limb. The 22mm female connector must be connected to the converter connected to the Y-piece.

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For insertion of MinimHal® on the inspiratory limb of an adult dual limb circuit:
The 15mm male connector must be connected to the converter itself connected to the inspiratory limb. The female connector must be connected to the Y-piece.

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For insertion of MinimHal® between the Y-piece and the endotracheal tube of an adult or pediatric dual limb circuit :
In case of an invasive ventilation, the male connector must be connected to the Y piece and the female connector must be connected to the converter. The converter must be connected to the patient's interface.
In case of a non-invasive ventilation, it is possible to connect the female connector directly to the patient's interface.

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
When using MinimHal® after the Y-piece, make sure that the dead volume (to which the chamber's dead volume is added) is compatible with the tidal volume delivered, and if necessary, disconnect MinimHal® from the circuit.

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DRUG ADMINISTRATION BY PRESSURIZED METERED DOSE INHALER

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Step 1
Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Remove the cap from the opening (c) and insert the canister (Fig.5).
Step 2
Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI (Fig.6).
Step 3
Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c) (Fig.7).

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The trademark Respimat® belongs to Boehringer Ingelheim. The trademark MinimHal® belongs to OptimHal-ProtectSom SAS.
Respimat® est une marque de Boehringer Ingelheim. MinimHal® est une marque d'OptimHal-ProtectSom SAS.
Brevets déposés conjointement par l'Université de Tours, l'INSERM et OptimHal-ProtectSom SAS.
Patents belongs to Tours University, INSERM, and OptimHal-ProtectSom SAS together.