
# CombiHaler® Solo

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
Inhalation chamber for mechanical ventilation circuit.

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


CombiHaler® Solo is an inhalation chamber intended for adult patient with invasive or non-invasive mechanical ventilation who have been prescribed aerosolized treatment from a pressurized metered-dose inhaler.

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- CombiHaler® Solo should be used by clinicians qualified in the care of mechanically ventilated patients.
- Please read the instructions for use of CombiHaler® Solo, the metered dose inhaler and keep them available at all times of use.

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- a. Male connector 22 mm
- b. Female connector 22 mm
- c. One opening for MDI
- d. One sealed cap
- e. One nozzle

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When putting the device into operation, it is imperative to ensure that:

- The device does not contain any foreign body, to avoid any risk of damage to the patient's airways or obstruction of the medical device.
- The cap is correctly positioned when MDI is not used.

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**Checking the operation of the nozzle**  
Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister (Fig.2). Operate the MDI and check that the nozzle correctly ensures its function (visually test with a puff of drugs) (Fig.3).

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**Insertion in the circuit**  
In a traditional adult adult dual limb of mechanical ventilation circuit, the CombiHaler® Solo must be inserted on the inspiratory limb.

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The arrow on the device must point towards the patient (Fig.4).

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If, for technical reason, CombiHaler® Solo would be inserted between the Y-piece and the patient, it's imperative that the inhalation chamber is removed from the system immediately after taking the medication and is never left in place while the circuit is being used.

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### DRUG ADMINISTRATION BY PRESSURIZED METERED DOSE INHALER

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**Step 1:**  
Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Remove the cap from the opening (c) and insert the canister (Fig.5).

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
**Step 2:**  
Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI (Fig.6).

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
**Step 3:**  
Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c) (Fig.7).

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## CLINICAL BENEFITS – DEVICE PERFORMANCE

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The chamber volume allows the propellant gas of the larger particles size to evaporate, producing smaller aerosol particles of the drug. Smaller particles improve the efficiency of aerosol delivery to the patient through the ventilation system. The clinical benefit is improved drug delivery.

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- The recommended shelf life of CombiHaler® Solo is that recommended for the mechanical ventilation system with which it is associated, not to exceed 7 days, or less if visual examination indicates a need for replacement.
- CombiHaler® Solo should be discarded along with the mechanical ventilation system in accordance with institutional guidelines.
- CombiHaler® Solo should never be used on more than one circuit and should be disposed of in the same manner as the circuit on which it was installed.
- Any total or partial replacement of the circuit must lead to the replacement of CombiHaler® Solo.
- CombiHaler® Solo must be stored in a dry place.

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- This device requires to be used by qualified staff to take care of patients who are under respiratory assistance.
- CombiHaler® Solo is a single patient product. Do not clean or sterilize. Do not reuse with other patients to avoid any risk of contamination.
- To ensure the holding of the device, CombiHaler® Solo must not be dislocated. In case of dislocation, CombiHaler® Solo must not be used, and must be replaced.
- CombiHaler® Solo is design to be inserted into an adult mechanical ventilation circuit.
- Any heat and moisture exchanger, located between the spacer and the patient, must be temporarily removed during the use of CombiHaler® Solo.
- Do not use CombiHaler® Solo with ventilation pressures exceeding 35 mbar, without permanent and careful attention.

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CombiHaler® Solo is a trademark of OptimHal-ProtectSom SAS.
CombiHaler® Solo est une marque de commerce d'OptimHal-ProtectSom SAS.
OptimHal-ProtectSom SAS